


Signia IF & Custom Hearing Aids – Replacing Batteries


Removing the battery

- Open the battery compartment.
- Remove the battery by hand.



Inserting the battery

- If the battery has a protective film, remove it only when you are ready to use the battery. To achieve maximum performance of the battery, once the film has been removed, refrain from putting the battery into the hearing aid for up to 2 minutes in order to allow it to activate.
- Insert the battery with the “+” symbol facing up.



- Carefully close the battery compartment. If you feel resistance, the battery is not inserted correctly.
- Do not attempt to close the battery compartment by force. It can be damaged.